



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05681 - Goat Anti-Tankyrase 2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TNKS2, tankyrase, TRF1-interacting ankyrin-related ADP-ribose polymerase 2, TNKL, TANK2, tankyrase 2, PARP-5b, PARP-5c, PARP5B, PARP5C

Official Symbol: TNKS2

Accession Number(s): NP_079511.1

Human GeneID(s): [80351](#)

Non-Human GeneID(s): 74493 (mouse)

Immunogen

Peptide with sequence SGRRCAGGGAACA, from the N Terminus of the protein sequence according to NP_079511.1.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Western Blot: Preliminary experiments gave no signal but low background in Human Muscle and Placenta extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow